# Tutoriel SnapGene Viewer

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Reverse complement

Create a sequence of 1000 N's

Create a placeholder file with no sequence

Topology:  
 Linear  Circular

File Name: Untitled 2.dna

Please enter a plain or formatted DNA sequence.

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Detect Common Features

Add detected features to Sequence\_test.dna:

Add	Feature	Location	Size	Type	Match
<input checked="" type="checkbox"/>	/gene	= lacZ fragment			
	/product	= LacZa fragment of β-galactosidase			
	/translation	= MTMITPLVRGLEIYRCMPWYPGARIRSFRRVDVIWIRYRRGGRS RTSGSIPNSPYSESYYNSLAVVLQRRDWENPGVTQLNRLAAHPPF ASWRNSEEARTDRPSQQLRSLNGEWRLMRYFLLTHLCGISHRIRQ SNHSTRPVAAH*			
		146 amino acids = 17,2 kDa			
<input checked="" type="checkbox"/>	<b>SK primer</b>	351 .. 367	17 bp	→ primer_bind	100.0%
	/note	= common sequencing primer, one of multiple similar variants			
<input checked="" type="checkbox"/>	<b>T7 promoter</b>	382 .. 400	19 bp	← promoter	100.0%
	/note	= promoter for bacteriophage T7 RNA polymerase			
<input checked="" type="checkbox"/>	<b>M13 fwd</b>	407 .. 423	17 bp	← primer_bind	100.0%
	/note	= common sequencing primer, one of multiple similar variants			
<input checked="" type="checkbox"/>	<b>ori</b>	2332 .. 2920	589 bp	→ rep_origin	99.8%
	/direction	= RIGHT			
	/note	= high-copy-number ColE1/pMB1/pBR322/pUC origin of replication			

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The screenshot displays the SnapGene Viewer interface for a circular DNA sequence of 3101 bp. The main window shows a circular map with various restriction sites and features. The sites are color-coded and labeled with their names and positions. Key features include the CAP binding site, lac operator, and several primers (M13 rev, SK primer, T7 promoter). The interface includes a menu bar (File, Edit, View, Enzymes, Features, Primers, Actions, Tools, Window, Help), a toolbar with icons for New, Open, Save, Print, Undo, Redo, Cut, Copy, and Paste, and a bottom panel with tabs for Map, Sequence, Enzymes, Features, Primers, and History. The Enzymes tab is currently selected, showing a list of 6+ unique cutters. The circular map is annotated with numerous restriction sites and features, including CAP binding site, lac operator, and various primers.

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The screenshot shows the SnapGene Viewer interface for a circular DNA sequence of 3101 bp. The main window displays the sequence in compact mode. The top menu bar includes File, Edit, View, Enzymes, Features, Primers, Actions, Tools, Window, and Help. Below the menu is a toolbar with icons for New, Open, Save, Print, Undo, Redo, Cut, Copy, and Paste. The sequence is shown in 5' to 3' orientation. The top strand is: `5' agcgccaatacgc... tcccgactggaagcgggcagtgagcgcaa`. The bottom strand is: `3' tcgcggttatg... tcaactcgcctattgta`. The sequence is annotated with several features: a CAP binding site (blue bar), a lac operator (blue bar), and a lacZa gene (purple bar). The amino acid translations are shown below the sequence. The top translation is: `Met Gln Leu Ala Arg Gln Val Ser Arg Leu Glu Ser Gly Gln`. The middle translation is: `Leu Ala Trp Tyr Ala Phe Arg Arg Glu Gly Arg Thr Pro Arg Asn Met Leu Ala Ala Pro Val Val Pro Lys Gly Val Pro Phe Arg Ala Thr Leu Ala Val`. The bottom translation is: `Cys Asn Ile His Thr Leu Glu Ser Met`. The bottom-most translation is: `Lys Val Cys Ser Val Ala Ile Val Met`. The interface also shows a Mini Map at the top right, a vertical toolbar on the left with options for salt, AC0TB, BTECA, CATAG, Asn, Arg, Ala, and CCA, and a status bar at the bottom with tabs for Map, Sequence, Enzymes, Features, Primers, and History.

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
The screenshot shows a DNA sequence in SnapGene Viewer. The top strand is 5' to 3' and the bottom strand is 3' to 5'. A blue selection box highlights the first 22 bases of the top strand. A purple box above the selection displays "219 .. 240 = 22 bp" and "Tm = 56-57°C". Below the sequence, a protein translation is shown with amino acid names: Met, Thr, Met, Ile, Thr, Pro, Leu, Val, Arg, Gly, Leu, Glu, Ile, Tyr, Arg, Cys. A purple arrow points to the start of the sequence.

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